

How to Collect Faeces Sample

※Be sure to read the instructions before collecting the faeces sample.

Precautions

- Do not tear the aluminum sheet.
- Be aware that if it is flushed together with a plenty of toilet paper, the toilet may become clogged.
- Be sure to take one faeces sample after one bowel movement.
- Be sure to take 2-day-samples at a short interval.
- Do not discard, drink, add water to the liquid in the sampling bottle.
- Do not obtain a faeces sample if you are menstruating.
- Do not use the sampling bottle directly on the body.
- Do not use the sampling bottle for any purpose other than obtaining a faeces sample.
- Be sure to keep it out of reach of children.
- Be careful to take the sample when use automatic cleaning function toilet.



Do not flush this paper into the toilet


Spread the collection paper inside the flush toilet and defecate onto the paper
 After obtaining the faeces sample, flush away the paper. However, be sure not to flush a plenty of papers at same time.



1 Complete details on label (Name, Gender, Age)




- ※Make sure to fill the collection date and time
- ※Make sure not to use Aqueous Pen



Twist the Green Cap on the sampling bottle to remove it

2

※Remember an excessive sample collecting cannot obtain an accurate examination



3

① Align the cap in sample direction as the sampling bottle


② Press it straight onto the bottle until a click sound is heard



※Do not turn in

Shut the cap until a Click sound is heard

※Make sure not to pull it again once pressed and shut the cap



※Store the samples at cool place

Place the two sampling bottles in submitting bag and submit immediately